
Gruppe

akkubetrieb

▶ = Info

⏪ = Eingabe

PEGA® TUBE (REF 10255)

**PEGA® TUBE**

CE 0123

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### Intended Use

PEGA® TUBE is a peristaltic tube, which is only suitable for use with PEGA® Infuse, PEGA® PCA, PEGA® plus and PEGA® light infusion pumps.

### Indication for Use

The PEGA® TUBE is suitable for intravenous, subcutaneous, epidural, and intrathecal administration of drugs from a PEGA® infusion pump.

View	No.	Description / Specification
	1	Connection (infusion tube) / Luer-Lock (male)
	2	Pressure tube (Silicone)
	3	Peristaltic tube (Silicone), Dead volume PEGA® TUBE $V_0 = 450 \mu\text{l}$
	4	Connection (reservoir) / Luer-Lock (male)

Fig. A: view PEGA® TUBE

### Conformity

Quality and safety of our products are the highest priority at Venner Medical. If a complaint should occur, please contact your responsible medical device consultant or

directly Venner Medical. In case of incidents with the PEGA products, report them to the competent authority in your country (e.g. BfArM in Germany).

## English

### Contraindication

Do not use drugs which are not compatible with the material (ABS, silicone rubber) of the PEGA® TUBE.

### Additional Information

PEGA® TUBE is latex, DEHP and pyrogen free.

### General Safety Instructions

- ◆ Please consult the manual of your PEGA® pump before using a PEGA® TUBE.
- ◆ This disposable item is intended to be used on a patient only once in order to guarantee optimum performance and hygiene.  
  
If the item is used again, even after reconditioning, any liability on the part of Venner Medical (Deutschland) GmbH shall cease to exist.
- ◆ Please also respect the manuals of the PEGA® BAG and PEGA® SPIKE.
- ◆ Carry out all of these steps on a clean and low germ concentration surface.
- ◆ Hands must be disinfected previously. Do not touch open luer-locks (danger of contamination).
- ◆ Use only original PEGA® accesso
- ◆ PEGA® TUBEs are provided in sterile packages. Do not use any disposable if the package appears damaged or opened.
- ◆ Use PEGA® TUBEs only once and discard the disposable as per local regulations.
- ◆ Only remove the vented caps of the inserted PEGA® TUBE when connecting to the PEGA® BAG and PEGA® COMBI LINE.
- ◆ PEGA® TUBE has to be substituted with every change of reservoir and must not be used longer than 30 days.
- ◆ A maximum delivery volume of 500 ml must not be exceeded.
- ◆ The pump with the PEGA® TUBE must not be positioned higher than 1 m above the heart of the patient during the therapy.
- ◆ The PEGA® pumps are correctly closed only if the clips are fully engaged. If the clips are not correctly engaged a safe operation of the infusion pumps cannot be warranted.

### Confirmation of sterilization

Yellow Color indicator for EO sterilization method on the transparent side of the sterile

pouch advises that the sterilization is done and passed.

**⚠ Warning!**

The product is designed for a maximum pressure of 1,800 hPa!


**Inserting**

Fig. B: view with infusion pump (example PEGA® Infuse)

1. Open the rear pump cover.
2. Open the locking clip.
3. Open the sterile packaging of the PEGA® TUBE.
4. Insert the peristaltic tube into the upper bearing of the pump. Ensure that the pressure tube is positioned inside the guides of the coloured pressure sensor.

5. Pull carefully on the Luer-lock-connection of the peristaltic tube so that this engages into the lower bearing of the pump.
6. If necessary, press the tube downwards with your finger so that it rests deep in the shaft.
7. Close the locking clip. The pressing block presses the peristaltic tube into the peristaltic drive.
8. Place the pump cover onto the pump and close it. The clips of the housing (one or two clips, depends on the PEGA® Pump version) must be heard to clip into place.

**⚠ Caution!**

If a PEGA® BAG is connected to the luer locks, make sure that the connections are always hand-tightened. This depends upon the position of the thread up to a three-quarter turn.

**Purchase Order**

Purchase only complete boxes.

REF.	Signification	Units per Box
10255	PEGA® TUBE	10

Table 1: purchase order

# English

## Meaning of the symbols on product and transport packaging

 <b>MD</b>	Medical Device	 <b>REF</b>	Article number
	Use by date		Open here
	Single use only. Do not re-use		Keep dry
	Do not re-sterilize		Temperature limit (+5° to +40°C)
	Do not use if package is damaged		Keep away from sunlight
	Caution or Warning		Does not contain phthalates (DEHP)
	Caution, consult instructions for use		Does not contain or no presence of natural rubber latex
	Date of manufacture and country of manufacture		Pyrogen-free
	Batch code		Manufacturer
	Sterilisation with ethylene oxide		
	Single sterile barrier system		
	CE mark with code of notified body - according to MDD 93/42 EEC and MDR		

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