Gebrauchsanweisung
Instructions for Use
Mode d'emploi

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umpe

kkubetrieb

► = Info

➔ = Eingabe

PEGA® SET 50 ml (REF 15055)
PEGA® SET 100 ml (REF 15105)
PEGA® SET 150 ml (REF 15155)
PEGA® SET 300 ml (REF 15305)

PEGA® SET

CE 0123

VENNER

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This page left intentionally blank

Intended Use

PEGA® SET is a complete, firmly connected catheter system consisting of a reservoir, a peristaltic tube and an infusion tube with a sterile filter for parenteral administration of drugs in ambulatory and in-patient infusion therapies.

PEGA® SET consists of a special check valve that prevents back-and freeflow.

PEGA® SET is only intended for use with PEGA® Infuse, PEGA® PCA, PEGA® plus and PEGA® light infusion pumps.

Indication for Use

The PEGA® SET is suitable for intravenous, subcutaneous, epidural, and intrathecal administration of drugs from a PEGA® infusion pump.


View	No.	Description / Specification
	1	Connection with siphon (patient) / Luer-Lock (male) Check valve
	2	0.2 µm Sterile filter with 0.1µm PVDF air-out-filter; with embossed device LOT number
	3	Infusion tube L = 1,050 mm, ID = 1 mm, Dead volume PEGA® LINE $V_0 = 1,400 \mu\text{l}$
	4	Segment of the pump (Silicone) consisting of
	5	Pressure tube (4) and Peristaltic tube (5)
	6	Connection (for filling) / Luer-Lock (female)
	7	Tube clamp (PVC)
	8	Medication reservoir 50 ml, 100 ml, 150 ml, 300 ml

Fig. A: view PEGA® SET

English

Contraindication

Do not use drugs which are not compatible with the material (PVC, PC, PP, PES, PVDF, PMMA, ABS, silicone rubber) of the PEGA® SET.

Additional Information

PEGA® SET is latex, DEHP and pyrogen free.

General Safety Instructions

- ◆ Please consult the manual of your PEGA® pump before using a PEGA® SET.
- ◆ This disposable item is intended to be used on a patient only once in order to guarantee optimum performance and hygiene.

If the item is used again, even after reconditioning, any liability on the part of Venner Medical (Deutschland) GmbH shall cease to exist.
- ◆ Carry out all of these steps on a clean and low germ concentration surface.
- ◆ Hands must be disinfected previously. Do not touch open luer-locks (danger of contamination).
- ◆ Use only original PEGA® accessories.
- ◆ PEGA® SETs are provided in sterile packages. Do not use any disposable if the package appears damaged or opened.

- ◆ Use PEGA® SETs only once and discard the disposable as per local regulations.
- ◆ Only open the vented caps of both luer-lock-connections on the reservoir when filling or priming the system.
- ◆ PEGA® SET must not be refilled and not be used longer than 30 days.
- ◆ The pump with the PEGA® SET must not be positioned higher than 1 m above the heart of the patient during the therapy.
- ◆ Do not exceed the indicated maximum volume for the PEGA® SET.
- ◆ If you want to store a filled PEGA® SET for a short time, please remove the vented cap on the luer-lock-connector on the patient's side and replace it with the cap provided.

Confirmation of sterilization

Yellow Color indicator for EO sterilization method on the transparent side of the sterile

pouch advises that the sterilization is done and passed.

⚠ Warning!

The product is designed for a maximum pressure of 1,800 hPa!

Filling the Reservoir


Fig. B: view PEGA® SET with syringe

1. Fill a luer- or a luer-lock-syringe with the medication solution.
2. Remove the filling cannula if necessary.
3. Remove the protective cap from the Luer-Lock connector on the reservoir filling tube and connect the syringe.
4. Inject the medication solution.
5. Prime the filled reservoir with the aid of the syringe.
6. Close the clamp on the reservoir filling tube.
7. Remove the syringe and fit the attached cap.

Preparation for use on the patient


Fig. C: View of pump (example PEGA® PCA) with PEGA® SET

1. Open the rear pump cover.
2. Open the locking clip.
3. Insert the peristaltic tube into the upper guide of the pump. Ensure that the pressure tube is positioned inside the guides of the coloured pressure sensor.
4. Pull carefully on the Luer-lock-connection of the peristaltic tube so that this engages into the lower guide of the pump.

English

5. If necessary, press the peristaltic tube down with your finger so that is positioned well into the guide.
6. Close the locking clip. The tube is fixed firmly in the peristaltic drive by the pressure block.
7. Place the pump lid on the pump and close it. The clips of the housing (one or two clips, depends on the PEGA® Pump version) must be heard to clip into place.

Caution

The pumps are correctly closed only if the clips are fully engaged. If the clips are not correctly engaged a safe operation of the PEGA® pumps cannot be warranted.

Priming

Before connecting it to the patient, prime the catheter system of the PEGA® SET following the instructions in the manual of your PEGA® pump.

Ensure that no air bubbles are left in the tubing or the filter on the patient side (ribbed

side). If this is not the case, flick the filter and tubing lightly with your finger until all the air bubbles are released.

Continue the priming procedure until a drop of medication has formed at the end of the catheter system.

Purchase Order

Purchase only complete boxes.

REF.	Signification	Units per Box
15055	PEGA® SET 50 ml	10
15105	PEGA® SET 100 ml	10
15155	PEGA® SET 150 ml	10
15305	PEGA® SET 300 ml	10

Table 1: purchase order

Conformity

Quality and safety of our products are the highest priority at Venner Medical. If a complaint should occur, please contact your responsible medical device consultant or

directly Venner Medical. In case of incidents with the PEGA products, report them to the competent authority in your country (e.g. BfArM in Germany).

Meaning of the symbols on product and transport packaging

 MD	Medical Device	 REF	Article number
	Use by date		Open here
	Single use only. Do not re-use		Keep dry
	Do not re-sterilize		Temperature limit (+5° to +40°C)
	Do not use if package is damaged		Keep away from sunlight
	Caution or Warning		Does not contain phthalates (DEHP)
	Caution, consult instructions for use		Does not contain or no presence of natural rubber latex
	Date of manufacture and country of manufacture		Pyrogen-free
 LOT	Batch code		Manufacturer
 STERILE EO	Sterilisation with ethylene oxide		
 SBS	Single sterile barrier system		
 CE 0123	CE mark with code of notified body - according to MDD 93/42 EEC and MDR		

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