
umpe

kkubetrieb

➔ = Info

➔ = Eingabe

PEGA® SPIKE (REF 10530)

PEGA® SPIKE

CE 0123

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This page left intentionally blank

Intended Use

PEGA® SPIKE is used as connection between PEGA® COMBI LINE and infusion bags with pierceable septum for parenteral administration of drugs in ambulatory and in-patient infusion therapies.

PEGA® SPIKE is only intended for use with the PEGA® Infuse, PEGA® PCA, PEGA® plus and PEGA® light infusion pumps.

PEGA® SPIKE comprises an unvented spike, an infusion tube and a luer lock connector with a safe lock (Luer-activated valve).

Indication for Use

The PEGA® SPIKE is suitable for intravenous, subcutaneous, epidural and intrathecal

administration of drugs from a PEGA® infusion pump.


View	No.	Description / Specification
	1	Connection (pump), Luer Lock (female), luer-activated valve, vented cap
	2	Infusion tube L = 300 mm, PVC ID = 2.2 mm, Dead volume PEGA® SPIKE: V_0 (with spike) = 1,700 μ l V_0 (without spike) = 1,400 μ l
	3	Spike (reservoir) (ABS), vented cap

Fig. A: view PEGA® SPIKE

English

Contraindication

Do not use drugs which are not compatible with the material (PVC, ABS, PC, silicone) of the PEGA® SPIKE.

Additional Information

PEGA® BAG is latex, DEHP and pyrogen free.

General Safety Instructions

- ◆ Please consult the manual of your PEGA® pump before using a PEGA® SPIKE.
- ◆ This disposable item is intended to be used on a patient only once in order to guarantee optimum performance and hygiene.

If the item is used again, even after reconditioning, any liability on the part of Venner Medical (Deutschland) GmbH shall cease to exist.
- ◆ Please also respect the manuals of the PEGA® COMBI LINE, PEGA® TUBE.
- ◆ Carry out all of these steps on a clean and low germ concentration surface.
- ◆ Hands must be disinfected previously. Do not touch open luer locks (danger of contamination).
- ◆ Use only original PEGA® accessories (with the exception of the infusion bag with pierceable septum).
- ◆ PEGA® SPIKEs are supplied in sterile packaging. Do not use if the package appears damaged or opened.
- ◆ Use PEGA® SPIKEs only once and dispose of in accordance with the disposal methods applicable in clinics and medical practices.
- ◆ PEGA® SPIKE must be replaced with every change of reservoir and must not be used for longer than 30 days.
- ◆ The maximum delivery volume of 500 ml must not be exceeded.
- ◆ The pump with the PEGA® SPIKE must not be positioned higher than 1 m above the heart of the patient during the therapy.
- ◆ Please ensure that the connection between the PEGA® SPIKE and the PEGA® COMBI LINE / PEGA® TUBE is sufficiently tight in order to avoid unintentional loosening of the connection.


Confirmation of sterilization

Yellow Color indicator for EO sterilization method on the transparent side of the sterile

pouch advises that the sterilization is done and passed.

⚠ Warning!

The product is designed for a maximum pressure of 1,800 hPa!

Connecting

1. Open the sterile packaging of the PEGA® SPIKE.
2. Remove the protective cap from the shaft of the PEGA® SPIKE.
3. Pierce the septum of the medication bag with the spike.
4. To remove the air, press carefully on the infusion bag until all the air has escaped through the PEGA® SPIKE.

Alternatively, you can also remove the air manually from the bag and the PEGA® SPIKE with the help of a syringe.

5. Remove the vented cap from the luer lock connector of the PEGA® SPIKE.
6. Connect the luer lock connections of the PEGA® SPIKE and PEGA® COMBI LINE carefully and screw the connection sufficiently tight.
7. Fit the complete system to the pump.

⚠ CAUTION!

Please refer to the instruction manuals for your PEGA® COMBI LINE.

Fig. B: PEGA® SPIKE and connection components

English

Priming

Before connecting it to the patient, prime the catheter system following the instructions in the manual of your PEGA® pump.

Ensure that no air bubbles are left in the tubing or the filter on the patient side (ribbed

side). If this is not the case, flick the filter and tubing lightly with your finger until all the air bubbles are released. Continue the priming procedure until a drop of medication has formed at the end of the catheter system.

Purchase Order

Purchase only complete boxes

REF.	Signification	Units per Box
10530	PEGA® SPIKE	10

Table 1: purchase order

Conformity

Quality and safety of our products are the highest priority at Venner Medical. If a complaint should occur, please contact your responsible medical device consultant or

directly Venner Medical. In case of incidents with the PEGA products, report them to the competent authority in your country (e.g. BfArM in Germany).

Meaning of the symbols on product and transport packaging

 MD	Medical Device	 REF	Article number
	Use by date		Open here
	Single use only. Do not re-use		Keep dry
	Do not re-sterilize		Temperature limit (+5° to +40°C)
	Do not use if package is damaged		Keep away from sunlight
	Caution or Warning		Does not contain phthalates (DEHP)
	Caution, consult instructions for use		Does not contain or no presence of natural rubber latex
	Date of manufacture and country of manufacture		Pyrogen-free
	Batch code		Manufacturer
	Sterilisation with ethylene oxide		
	Single sterile barrier system		
	CE mark with code of notified body - according to MDD 93/42 EEC and MDR		

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