


# PEGA<sup>®</sup> INFUSE

SMALL. SIMPLE. FLEXIBLE.




# Simply, simple.



With only 3 buttons the PEGA® Infuse is very easy to use. The intuitive design enables you to operate the infusion pump safely, supporting your daily work and reducing the need for training, allowing you more time for other tasks.



Simply,  
flexible.



**Customise your PEGA® Infuse!**

In addition to individual therapy settings, the PEGA® Infuse allows you to quickly and easily pre-set, name and run up to six drug protocols. This reduces programming errors and provides additional safety for users and patients.




Simply,  
automated.



**Designed for your modern pain therapy!**

The automatic bolus is administered at defined intervals. Simply set the desired interval time and start your therapy. This reduces the time and effort required for pain therapy, leading to a quick mobilization and a satisfied patient.






**Even more time!**

Simply charge the PEGA® Infuse during therapy. Connect the Venner approved Powerbank to the PEGA® Infuse. Thus, extending the running time of the pump, no matter where you are.



Be mobile,  
stay mobile.





**Small & lightweight**

The ideal companion for the mobile patient.

**Modern design**


Proven quality in a new design.

**Full range of disposables**

Suitable for various infusion therapy protocols.

**Integrated battery**


To protect the environment the PEGA® Infuse eliminates the need for disposable batteries.







PEGA Infuse






# Simply, more.

The intuitive, advanced software allows user-friendly and code-protected programming of the PEGA® Infuse while offering you multiple therapy options. Choose between therapies for continuous infusion, PCA infusion, intermittent bolus or clinician activated bolus.

Use the PEGA® Infuse PumpReader software, exclusively developed for the PEGA® Infuse, to read out and document your PEGA® Infuse therapies.





# PEGA<sup>®</sup> Infuse

**Specially designed with over 30 years of experience in infusion technology.**

With more than 30 years of experience Venner Medical develops and distributes small, easy to use and reliable infusion pumps for stationary and ambulatory infusion therapy.

We want our products to improve the quality of life of all people. This is why we put great emphasis on the safety and reliability of our products.

We work closely with our customers and patients to develop innovative products and future therapy concepts.

For the best possible support in the hospital and at home, our specialised teams are there to assist you. Whether it is technical questions, personal service or customised training and instruction - we are here for you.



Anaesthesia.  
Intensive care.  
Emergency Medicine.  
Orthopaedics.  
Obstetrics.



Palliative care.  
Pain therapy.  
Chemotherapy.  
Desferal therapy.





## Suitable accessories:

CE

<b>Product Code</b>	<b>Description</b>
<b>60100</b>	PEGA® Infuse BOX 100 (≤ 100 ml)
<b>60150</b>	PEGA® Infuse BOX 150 (≤ 150 ml)
<b>60810</b>	PEGA® Powerbank
<b>60910</b>	PEGA® POUCH S Single-Use Pouch
<b>60930</b>	PEGA® POUCH L Single-Use Pouch

## Suitable accessories:

<b>Product Code</b>	<b>Description</b>
<b>60180</b>	PEGA® Infuse Box Key
<b>60104</b>	PEGA® Infuse PumpReader

# Order details

PEGA® Infuse Carry Case Content:

CE 0123

Amount	Content
1	PEGA® Infuse
1	PEGA® Infuse BOLUS BUTTON (REF:60211)
1	PEGA® POWER SUPPLY USB (REF:60830)
1	Instruction for Use
1	Short Instruction PEGA® Infuse (Doctor)
1	Short Instruction PEGA® Infuse (Healthcare)

# Order details

PEGA® Consumables:

CE 0123

Product Code	Description	Qty./box
14150	PEGA® BAG 150 ml	10
14300	PEGA® BAG 300 ml	10
15055	PEGA® SET 50 ml	10
15105	PEGA® SET 100 ml	10
15155	PEGA® SET 150 ml	10
15305	PEGA® SET 300 ml	10
18050	PEGA® SET curled 50 ml	10
18100	PEGA® SET curled 100 ml	10
18150	PEGA® SET curled 150 ml	10

CE 0123

<b>Product code</b>	<b>Description</b>	<b>Qty./box</b>
<b>18300</b>	PEGA® SET curled 300 ml	<b>10</b>
<b>10521S</b>	PEGA® COMBI LINE S	<b>10</b>
<b>10522E</b>	PEGA® COMBI LINE E	<b>10</b>
<b>10526E</b>	PEGA® COMBI LINE E with NRFit® connector	<b>10</b>
<b>10523Y</b>	PEGA® COMBI LINE Y	<b>10</b>
<b>10530</b>	PEGA® SPIKE	<b>10</b>
<b>14050</b>	PEGA® BAG 50 ml	<b>10</b>
<b>14100</b>	PEGA® BAG 100 ml	<b>10</b>



 **Manufacturer:**

Venner Medical (Deutschland) GmbH

Gut Uhlenhorst 4

24229 Dänischenhagen

Phone: +49 (0) 4349 - 9154 - 0

Fax: +49 (0) 4349 - 9154 - 10

Email: [info@vennermedical.de](mailto:info@vennermedical.de)

[www.vennermedical.de/en](http://www.vennermedical.de/en)



Find out more